

Compound TN12034

Chemical Properties


CAS No. :

Formula: C32H38O19

Molecular Weight: 726.64

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

methyl (2S,3Z,4S)-3-ethylidene-4-(2-oxo-2-(((2R,3S,4S,5R,6S)-3,4,5-trihydroxy-6-((6-hydroxy-2-oxo-2H-chromen-7-yl)oxy)oxan-2-yl]methoxy)ethyl)-2-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)-3,4-dihydro-2H-pyran-5-carboxylate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3762 mL	6.881 mL	13.762 mL
5 mM	0.2752 mL	1.3762 mL	2.7524 mL
10 mM	0.1376 mL	0.6881 mL	1.3762 mL
50 mM	0.0275 mL	0.1376 mL	0.2752 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481