

## Compound TN12072

## Chemical Properties


CAS No. :

Formula: C31H40O15

Molecular Weight: 652.65

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description

5-[[3,4-dihydroxy-4-(hydroxymethyl)oxolan-2-yl]oxy]-2-(hydroxymethyl)-6-[4-[3-(hydroxymethyl)-5-[(1E)-3-hydroxyprop-1-en-1-yl]-7-methoxy-2,3-dihydro-1-benzofuran-2-yl]-2-methoxyphenoxy]oxane-3,4-diol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5322 mL	7.6611 mL	15.3221 mL
5 mM	0.3064 mL	1.5322 mL	3.0644 mL
10 mM	0.1532 mL	0.7661 mL	1.5322 mL
50 mM	0.0306 mL	0.1532 mL	0.3064 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481