

## 3-[(3-[(2R,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy)decanoyl]oxy]decanoic acid

### Chemical Properties


CAS No. :

Formula: C<sub>26</sub>H<sub>48</sub>O<sub>9</sub>

Molecular Weight: 504.66

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description	3-[(3-[(2R,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy)decanoyl]oxy]decanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9815 mL	9.9077 mL	19.8153 mL
5 mM	0.3963 mL	1.9815 mL	3.9631 mL
10 mM	0.1982 mL	0.9908 mL	1.9815 mL
50 mM	0.0396 mL	0.1982 mL	0.3963 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481