

2-(3,4-dihydroxyphenyl)-3,3,5,7-tetrahydroxy-2-methoxy-3,4-dihydro-2H-1-benzopyran-4-one

Chemical Properties


CAS No. :

Formula: C₁₆H₁₄O₉

Molecular Weight: 350.28

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-(3,4-dihydroxyphenyl)-3,3,5,7-tetrahydroxy-2-methoxy-3,4-dihydro-2H-1-benzopyran-4-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8549 mL	14.2743 mL	28.5486 mL
5 mM	0.571 mL	2.8549 mL	5.7097 mL
10 mM	0.2855 mL	1.4274 mL	2.8549 mL
50 mM	0.0571 mL	0.2855 mL	0.571 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481