

Compound TN12138

Chemical Properties


CAS No. :

Formula: C₂₃H₂₄O₁₃

Molecular Weight: 508.43

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	[(2R,3S,4S,5R,6R)-6-[(5,8-dihydroxy-4-oxo-1,2,3,4-tetrahydronaphthalen-1-yl)oxy]-3,4,5-trihydroxyoxan-2-yl]methyl 3,4,5-trihydroxybenzoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9668 mL	9.8342 mL	19.6684 mL
5 mM	0.3934 mL	1.9668 mL	3.9337 mL
10 mM	0.1967 mL	0.9834 mL	1.9668 mL
50 mM	0.0393 mL	0.1967 mL	0.3934 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481