

Compound TN12157

Chemical Properties


CAS No. :

Formula: C₁₈H₂₈O₁₁

Molecular Weight: 420.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

methyl (1S,4aS,8R)-6-methoxy-8-methyl-1-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl)oxy)-1H,4aH,5H,6H,8H,8aH-pyrano[3,4-c]pyran-4-carboxylate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3786 mL	11.8932 mL	23.7863 mL
5 mM	0.4757 mL	2.3786 mL	4.7573 mL
10 mM	0.2379 mL	1.1893 mL	2.3786 mL
50 mM	0.0476 mL	0.2379 mL	0.4757 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481