

Compound TN12175

Chemical Properties


CAS No. :

Formula: C₂₉H₄₆O₈

Molecular Weight: 522.68

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-[3-acetyl-6,9a,9b-trimethyl-7-(prop-1-en-2-yl)-4-[[[(2S,3R,4S,5S)-3,4,5-trihydroxyoxan-2-yl]oxy]-dodecahydro-1H-cyclopenta[a]naphthalen-6-yl]propanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9132 mL	9.5661 mL	19.1322 mL
5 mM	0.3826 mL	1.9132 mL	3.8264 mL
10 mM	0.1913 mL	0.9566 mL	1.9132 mL
50 mM	0.0383 mL	0.1913 mL	0.3826 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481