

## Compound TN12177

## Chemical Properties


CAS No. :

Formula: C<sub>36</sub>H<sub>60</sub>O<sub>9</sub>

Molecular Weight: 636.87

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-{3-[(4E)-2-hydroxy-6-methoxy-6-methylhept-4-en-2-yl]-6,9a,9b-trimethyl-7-(prop-1-en-2-yl)-4-[[[(2S,3R,4S,5S)-3,4,5-trihydroxyoxan-2-yl]oxy]-dodecahydro-1H-cyclopenta[a]naphthalen-6-yl}propanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5702 mL	7.8509 mL	15.7018 mL
5 mM	0.314 mL	1.5702 mL	3.1404 mL
10 mM	0.157 mL	0.7851 mL	1.5702 mL
50 mM	0.0314 mL	0.157 mL	0.314 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481