

## 2-hydroxy-4-methoxy-6-[(1E)-2-phenylethenyl]benzoic acid

## Chemical Properties


CAS No. :

Formula: C<sub>16</sub>H<sub>14</sub>O<sub>4</sub>

Molecular Weight: 270.28

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description

2-hydroxy-4-methoxy-6-[(1E)-2-phenylethenyl]benzoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6999 mL	18.4993 mL	36.9987 mL
5 mM	0.740 mL	3.6999 mL	7.3997 mL
10 mM	0.370 mL	1.8499 mL	3.6999 mL
50 mM	0.074 mL	0.370 mL	0.740 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481