

## Compound TN12194

## Chemical Properties


CAS No. :

Formula: C<sub>20</sub>H<sub>22</sub>O<sub>9</sub>

Molecular Weight: 406.39

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(2R,3R)-7-(((2S,3R,4R,5S)-3,4-dihydroxy-5-(hydroxymethyl)oxolan-2-yl]oxy)-2-(4-hydroxyphenyl)-3,4-dihydro-2H-1-benzopyran-3,5-diol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4607 mL	12.3035 mL	24.6069 mL
5 mM	0.4921 mL	2.4607 mL	4.9214 mL
10 mM	0.2461 mL	1.2303 mL	2.4607 mL
50 mM	0.0492 mL	0.2461 mL	0.4921 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481