

Compound TN12196

Chemical Properties


CAS No. :

Formula: C₄₁H₇₀O₁₂

Molecular Weight: 755

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|-------------|--|
| Description | (2R,3R,4S,5S,6R)-2-(((3aR,3bR,5aR,7S,9aR,11R,11aR)-7-(((2R,3R,4R,5S)-3,4-dihydroxy-5-(hydroxymethyl)oxolan-2-yl]oxy)-1-[(2S,5R)-5-(2-hydroxypropan-2-yl)-2-methyloxolan-2-yl]-3a,3b,6,6,9a-pentamethyl-hexadecahydro-1H-cyclopenta[a]phenanthren-11-yl]oxy)-6-methyloxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research. |
|-------------|--|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|-----------|
| 1 mM | 1.3245 mL | 6.6225 mL | 13.245 mL |
| 5 mM | 0.2649 mL | 1.3245 mL | 2.649 mL |
| 10 mM | 0.1325 mL | 0.6623 mL | 1.3245 mL |
| 50 mM | 0.0265 mL | 0.1325 mL | 0.2649 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481