

Compound TN12215

Chemical Properties


CAS No. :

Formula: C40H68O12

Molecular Weight: 740.97

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-((1-[5-(2-hydroxypropan-2-yl)-2-methyloxolan-2-yl]-3a,3b,6,6,9a-pentamethyl-11-[(3,4,5-trihydroxyoxan-2-yl)oxy]-hexadecahydro-1H-cyclopenta[a]phenanthren-7-yl)oxy)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3496 mL	6.7479 mL	13.4958 mL
5 mM	0.2699 mL	1.3496 mL	2.6992 mL
10 mM	0.135 mL	0.6748 mL	1.3496 mL
50 mM	0.027 mL	0.135 mL	0.2699 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481