


2-{2-[2-(2-acetamidopropanamido)-3-hydroxybutanamido]-3-methylbutanamido}pentanedioic acid

Chemical Properties

CAS No. :

Formula: C19H32N4O9

Molecular Weight: 460.48



Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

2-{2-[2-(2-acetamidopropanamido)-3-hydroxybutanamido]-3-methylbutanamido}pentanedioic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1716 mL	10.8582 mL	21.7165 mL
5 mM	0.4343 mL	2.1716 mL	4.3433 mL
10 mM	0.2172 mL	1.0858 mL	2.1716 mL
50 mM	0.0434 mL	0.2172 mL	0.4343 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481