

## Compound TN12227

## Chemical Properties


CAS No. :

Formula: C53H90O23

Molecular Weight: 1095.28

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

## Description

(2S,3R,4S,5S,6R)-2-[[[(2S)-2-[(1S,3aR,3bR,7R,8R,9aR,11R)-7-[[[(2R,3R,4S,5S,6R)-4,5-dihydroxy-3-[[[(2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]-6-[[[(3,4,5-trihydroxyoxan-2-yl)oxy]methyl]oxan-2-yl]oxy]-8,11-dihydroxy-3a,3b,6,6,9a-pentamethyl-hexadecahydro-1H-cyclopenta[a]phenanthren-1-yl]-6-methylhept-5-en-2-yl]oxy]-6-(hydroxymethyl)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.913 mL	4.565 mL	9.1301 mL
5 mM	0.1826 mL	0.913 mL	1.826 mL
10 mM	0.0913 mL	0.4565 mL	0.913 mL
50 mM	0.0183 mL	0.0913 mL	0.1826 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481