

Compound TN12229

Chemical Properties


CAS No. :

Formula: C31H52O5

Molecular Weight: 504.75

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	methyl 3-{3-[(4E)-2,6-dihydroxy-6-methylhept-4-en-2-yl]-4-hydroxy-6,9a,9b-trimethyl-7-(prop-1-en-2-yl)-dodecahydro-1H-cyclopenta[a]naphthalen-6-yl}propanoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9812 mL	9.9059 mL	19.8118 mL
5 mM	0.3962 mL	1.9812 mL	3.9624 mL
10 mM	0.1981 mL	0.9906 mL	1.9812 mL
50 mM	0.0396 mL	0.1981 mL	0.3962 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481