

Compound TN12240

Chemical Properties


CAS No. :

Formula: C53H90O24

Molecular Weight: 1111.28

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

(2S,3R,4S,5S,6R)-2-[[[(2S,4E)-2-[(1S,3aR,3bR,7R,8R,9aR,11R)-7-[[[(2R,3R,4S,5S,6R)-4,5-dihydroxy-6-(hydroxymethyl)-3-[[[(2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]oxan-2-yl]oxy]-8,11-dihydroxy-3a,3b,6,6,9a-pentamethyl-hexadecahydro-1H-cyclopenta[a]phenanthren-1-yl]-6-hydroxy-6-methylhept-4-en-2-yl]oxy]-6-[[[(2S,3R,4S,5R)-3,4,5-trihydroxyoxan-2-yl]oxy)methyl]oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8999 mL	4.4993 mL	8.9986 mL
5 mM	0.180 mL	0.8999 mL	1.7997 mL
10 mM	0.090 mL	0.4499 mL	0.8999 mL
50 mM	0.018 mL	0.090 mL	0.180 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481