

Compound TN12243

Chemical Properties


CAS No. :

Formula: C₃₆H₆₀O₁₀

Molecular Weight: 652.87

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	methyl 3-[(4R,6S,7S,9aR,9bR)-3-[(2S,4E)-2-[[[(2S,3R,4R,5S)-3,4-dihydroxy-5-(hydroxymethyl)oxolan-2-yl]oxy]-6-hydroperoxy-6-methylhept-4-en-2-yl]-4-hydroxy-6,9a,9b-trimethyl-7-(prop-1-en-2-yl)-dodecahydro-1H-cyclopenta[a]naphthalen-6-yl]propanoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5317 mL	7.6585 mL	15.317 mL
5 mM	0.3063 mL	1.5317 mL	3.0634 mL
10 mM	0.1532 mL	0.7658 mL	1.5317 mL
50 mM	0.0306 mL	0.1532 mL	0.3063 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481