

Compound TN12258

Chemical Properties


CAS No. :

Formula: C₂₆H₃₀O₆

Molecular Weight: 438.52

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

3-[5-hydroxy-7-methoxy-2,3-dimethyl-6-(3-methylbut-2-en-1-yl)-4-oxo-3,4-dihydro-2H-1-benzopyran-8-yl]-3-phenylpropanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2804 mL	11.402 mL	22.804 mL
5 mM	0.4561 mL	2.2804 mL	4.5608 mL
10 mM	0.228 mL	1.1402 mL	2.2804 mL
50 mM	0.0456 mL	0.228 mL	0.4561 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481