

Compound TN12263

Chemical Properties


CAS No. :

Formula: C₂₁H₃₂O₁₂

Molecular Weight: 476.48

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2S,3R,4R,5R,6S)-2-(((2R,3R,4S,5R,6R)-3,5-dihydroxy-2-[2-(4-hydroxy-3-methoxyphenyl)ethoxy]-6-(hydroxymethyl)oxan-4-yl]oxy)-6-methyloxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0987 mL	10.4936 mL	20.9872 mL
5 mM	0.4197 mL	2.0987 mL	4.1974 mL
10 mM	0.2099 mL	1.0494 mL	2.0987 mL
50 mM	0.042 mL	0.2099 mL	0.4197 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481