

Compound TN12266

Chemical Properties


CAS No. :

Formula: C37H60O9

Molecular Weight: 648.88

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-{3-[2-((5-[(acetyloxy)methyl]-3,4-dihydroxyoxolan-2-yl)oxy)-6-methylhept-5-en-2-yl]-4-hydroxy-6,9a,9b-trimethyl-7-(prop-1-en-2-yl)-dodecahydro-1H-cyclopenta[a]naphthalen-6-yl}propanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5411 mL	7.7056 mL	15.4112 mL
5 mM	0.3082 mL	1.5411 mL	3.0822 mL
10 mM	0.1541 mL	0.7706 mL	1.5411 mL
50 mM	0.0308 mL	0.1541 mL	0.3082 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481