

5-(4-hydroxyphenyl)-3-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]pentanoic acid

Chemical Properties


CAS No. :

Formula: C₁₇H₂₄O₉

Molecular Weight: 372.37

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-(4-hydroxyphenyl)-3-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]pentanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6855 mL	13.4275 mL	26.855 mL
5 mM	0.5371 mL	2.6855 mL	5.371 mL
10 mM	0.2686 mL	1.3428 mL	2.6855 mL
50 mM	0.0537 mL	0.2686 mL	0.5371 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481