

## Compound TN12335

## Chemical Properties


CAS No. :

Formula: C58H90O22

Molecular Weight: 1139.34

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

## Description

(1S,3aS,7S,9aS,9bS,10S,11S,11aS)-1-acetyl-10-(acetyloxy)-9a-[(acetyloxy)methyl]-7-[[[(2R,4S,5R,6R)-5-[[[(2S,4S,5R,6R)-5-[[[(2S,4S,5R,6R)-5-[[[(2S,3R,4S,5R,6R)-3,5-dihydroxy-4-methoxy-6-methyloxan-2-yl]oxy]-4-methoxy-6-methyloxan-2-yl]oxy]-4-methoxy-6-methyloxan-2-yl]oxy]-4-methoxy-6-methyloxan-2-yl]oxy]-3a-hydroxy-11a-methyl-1H,2H,3H,3aH,3bH,4H,6H,7H,8H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-11-yl (2E)-2-methylbut-2-enoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8777 mL	4.3885 mL	8.777 mL
5 mM	0.1755 mL	0.8777 mL	1.7554 mL
10 mM	0.0878 mL	0.4389 mL	0.8777 mL
50 mM	0.0176 mL	0.0878 mL	0.1755 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481