

## Compound TN12341

## Chemical Properties


CAS No. :

Formula: C48H80O19

Molecular Weight: 961.15

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(3aR,3bR,9aR)-9-hydroxy-3a,3b,6,6,9a-pentamethyl-1-[(2R)-6-methyl-2-[[[(2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]hept-5-en-2-yl]-7-[[[(2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-[[[(2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]methyl]oxan-2-yl]oxy]-hexadecahydro-1H-cyclopenta[a]phenanthren-11-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0404 mL	5.2021 mL	10.4042 mL
5 mM	0.2081 mL	1.0404 mL	2.0808 mL
10 mM	0.104 mL	0.5202 mL	1.0404 mL
50 mM	0.0208 mL	0.104 mL	0.2081 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481