

Compound TN12350

Chemical Properties


CAS No. :

Formula: C₁₉H₃₀O₈

Molecular Weight: 386.44

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

(1R,5S,8R)-1,5-dimethyl-8-[(1E,3R)-3-[(2S,3S,4R,5R,6S)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]but-1-en-1-yl]-6-oxabicyclo[3.2.1]octan-3-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5877 mL	12.9386 mL	25.8772 mL
5 mM	0.5175 mL	2.5877 mL	5.1754 mL
10 mM	0.2588 mL	1.2939 mL	2.5877 mL
50 mM	0.0518 mL	0.2588 mL	0.5175 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481