

2-(4-hydroxy-3-methoxyphenyl)-4-[(hydroxy-3-methoxyphenyl)methyl]-3-(hydroxymethyl)oxolane-3,4-diol


## Chemical Properties

CAS No. :

Formula: C<sub>20</sub>H<sub>24</sub>O<sub>8</sub>

Molecular Weight: 392.4

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	2-(4-hydroxy-3-methoxyphenyl)-4-[(4-hydroxy-3-methoxyphenyl)methyl]-3-(hydroxymethyl)oxolane-3,4-diol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5484 mL	12.7421 mL	25.4842 mL
5 mM	0.5097 mL	2.5484 mL	5.0968 mL
10 mM	0.2548 mL	1.2742 mL	2.5484 mL
50 mM	0.051 mL	0.2548 mL	0.5097 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481