

methyl 2-[1,5-dihydroxy-4a,8-dimethyl-1,2,3,4,4a,5,6,8a-octahydronaphthalen-2-yl]prop-2-enoate


Chemical Properties

CAS No. :

Formula: C₁₆H₂₄O₄

Molecular Weight: 280.36

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	methyl 2-[(1S,2S,4aR)-1,5-dihydroxy-4a,8-dimethyl-1,2,3,4,4a,5,6,8a-octahydronaphthalen-2-yl]prop-2-enoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5668 mL	17.8342 mL	35.6684 mL
5 mM	0.7134 mL	3.5668 mL	7.1337 mL
10 mM	0.3567 mL	1.7834 mL	3.5668 mL
50 mM	0.0713 mL	0.3567 mL	0.7134 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481