

## Compound TN12359

## Chemical Properties


CAS No. :

Formula: C58H94O25

Molecular Weight: 1191.36

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(2R,3S,4S,5R,6S)-6-(((2R,3R,4S,6R)-6-(((7S,9aR,11aS)-1-((1R)-1-(((2R,3R,4S,5S,6R)-5-hydroxy-4-methoxy-6-methyl-3-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(((2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)methyl)oxan-2-yl]oxy)oxan-2-yl]oxy)ethyl]-9a,11a-dimethyl-1H,2H,3H,3aH,3bH,4H,6H,7H,8H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-7-yl]oxy)-4-methoxy-2-methyloxan-3-yl]oxy)-5-(acetyloxy)-4-methoxy-2-methyloxan-3-yl acetate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8394 mL	4.1969 mL	8.3938 mL
5 mM	0.1679 mL	0.8394 mL	1.6788 mL
10 mM	0.0839 mL	0.4197 mL	0.8394 mL
50 mM	0.0168 mL	0.0839 mL	0.1679 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481