

## Compound TN12362

## Chemical Properties


CAS No. :

Formula: C<sub>51</sub>H<sub>80</sub>O<sub>23</sub>

Molecular Weight: 1061.18

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

## Description

(1R,3aR,5aS,7S,9aS,11aR)-3a-hydroxy-7-(((2R,4S,5R,6R)-5-(((2S,3R,4R,5R,6R)-3-hydroxy-4-methoxy-6-methyl-5-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(((2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)methyl)oxan-2-yl]oxy)oxan-2-yl]oxy)-4-methoxy-6-methyl-oxan-2-yl]oxy)-9a,11a-dimethyl-1-(5-oxo-2,5-dihydrofuran-3-yl)-hexadecahydro-1H-cyclopenta[a]phenanthren-2-yl acetate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9423 mL	4.7117 mL	9.4235 mL
5 mM	0.1885 mL	0.9423 mL	1.8847 mL
10 mM	0.0942 mL	0.4712 mL	0.9423 mL
50 mM	0.0188 mL	0.0942 mL	0.1885 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481