

## Compound TN12393

## Chemical Properties


CAS No. :

Formula: C<sub>64</sub>H<sub>104</sub>O<sub>30</sub>

Molecular Weight: 1353.51

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

## Description

(2R,3S,4S,5R,6S)-5-(acetyloxy)-6-(((2R,3R,4S,6R)-6-(((7S)-1-((1R)-1-(((2R,3R,4S,5S,6R)-5-hydroxy-4-methoxy-6-methyl-3-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(((2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-((hydroxymethyl)oxan-2-yl]oxy)methyl)oxan-2-yl]oxy)methyl)oxan-2-yl]oxy)oxan-2-yl]oxy)ethyl]-9a,11a-dimethyl-1H,2H,3H,3aH,3bH,4H,6H,7H,8H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-7-yl]oxy)-4-methoxy-2-methyloxan-3-yl]oxy)-4-methoxy-2-methyloxan-3-yl acetate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7388 mL	3.6941 mL	7.3882 mL
5 mM	0.1478 mL	0.7388 mL	1.4776 mL
10 mM	0.0739 mL	0.3694 mL	0.7388 mL
50 mM	0.0148 mL	0.0739 mL	0.1478 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481