

## Compound TN12407

## Chemical Properties


CAS No. :

Formula: C<sub>65</sub>H<sub>106</sub>O<sub>31</sub>

Molecular Weight: 1383.53

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

## Description

(2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(((2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)methyl)oxan-2-yl (4aS,6aS,6bR,10S,12aR,14bS)-10-(((2S,3R,4S,5S)-3-(((2S,3R,4R,5S,6S)-4-(((2S,3R,4R,5S,6R)-3,4-dihydroxy-6-(hydroxymethyl)-5-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)oxan-2-yl]oxy)-3,5-dihydroxy-6-methyloxan-2-yl]oxy)-4,5-dihydroxyoxan-2-yl]oxy)-2,2,6a,6b,9,9,12a-heptamethyl-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,12b,13,14b-icosahydropicene-4a-carboxylate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7228 mL	3.6139 mL	7.2279 mL
5 mM	0.1446 mL	0.7228 mL	1.4456 mL
10 mM	0.0723 mL	0.3614 mL	0.7228 mL
50 mM	0.0145 mL	0.0723 mL	0.1446 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481