

Compound TN12431

Chemical Properties


CAS No. :

Formula: C27H34O14

Molecular Weight: 582.56

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2,3-dihydroxy-1-(3-methoxy-4-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]phenyl)propyl (2E)-3-(4-hydroxy-3,5-dimethoxyphenyl)prop-2-enoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7166 mL	8.5828 mL	17.1656 mL
5 mM	0.3433 mL	1.7166 mL	3.4331 mL
10 mM	0.1717 mL	0.8583 mL	1.7166 mL
50 mM	0.0343 mL	0.1717 mL	0.3433 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481