

## Compound TN12432

## Chemical Properties


CAS No. :

Formula: C<sub>40</sub>H<sub>42</sub>O<sub>12</sub>

Molecular Weight: 714.76

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2,3-bis[(4-hydroxy-3-methoxyphenyl)methyl]-4-(((2E)-3-(4-hydroxy-3-methoxyphenyl)prop-2-enyl)oxy)butyl (2E)-3-(4-hydroxy-3-methoxyphenyl)prop-2-enoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3991 mL	6.9954 mL	13.9907 mL
5 mM	0.2798 mL	1.3991 mL	2.7981 mL
10 mM	0.1399 mL	0.6995 mL	1.3991 mL
50 mM	0.028 mL	0.1399 mL	0.2798 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481