

## Compound TN12450

## Chemical Properties


CAS No. :

Formula: C61H96O33

Molecular Weight: 1357.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

## Description

(2S,3R,4S,5R,6R)-2-[(1'S,2S,3S,4S,4'S,7'R,9'S,13'R,14'R,16'R)-14'-[[[(2R,3R,4S,5R,6R)-3-[[[(2S,3R,4R,5S,6S)-4-[[[(2S,3R,4R)-3,4-dihydroxy-4-(hydroxymethyl)oxolan-2-yl]oxy]-3,5-dihydroxy-6-methyloxan-2-yl]oxy]-5-hydroxy-6-(hydroxymethyl)-4-[[[(2S,3R,4S,5R)-3,4,5-trihydroxyoxan-2-yl]oxy]oxan-2-yl]oxy]-9',13'-dimethyl-5-methylidene-7'-[[[(2R,3R,4S,5R,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]methyl]-5'-oxaspiro[oxane-2,6'-pentacyclo[10.8.0.0<sup>2,9</sup>].0<sup>4,8</sup>].0<sup>13,18</sup>]]icosan]-18'-ene-3,16'-dioloxy]-6-methyloxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7367 mL	3.6835 mL	7.367 mL
5 mM	0.1473 mL	0.7367 mL	1.4734 mL
10 mM	0.0737 mL	0.3683 mL	0.7367 mL
50 mM	0.0147 mL	0.0737 mL	0.1473 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481