

Compound TN12453

Chemical Properties


CAS No. :

Formula: C61H96O32

Molecular Weight: 1341.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2S,3R,4S,5R,6R)-2-[(1'S,2S,3S,4S,4'S,7'R,9'S,13'R,14'R,16'R)-14'-[[[(2R,3R,4S,5R,6R)-3-[[[(2S,3R,4R,5S,6S)-4-[[[(2S,3R,4R)-3,4-dihydroxy-4-(hydroxymethyl)oxolan-2-yl]oxy]-3,5-dihydroxy-6-methyloxan-2-yl]oxy]-5-hydroxy-6-methyl-4-[[[(2S,3R,4S,5R)-3,4,5-trihydroxyoxan-2-yl]oxy]oxan-2-yl]oxy]-9',13'-dimethyl-5-methylidene-7'-[[[(2R,3R,4S,5R,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]methyl]-5'-oxaspiro[oxane-2,6'-pentacyclo[10.8.0.0 ^{2,9}].0 ^{4,8}].0 ^{13,18}]]icosan]-18'-ene-3,16'-dioloxyl]-6-methyloxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7455 mL	3.7274 mL	7.4548 mL
5 mM	0.1491 mL	0.7455 mL	1.491 mL
10 mM	0.0745 mL	0.3727 mL	0.7455 mL
50 mM	0.0149 mL	0.0745 mL	0.1491 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481