

5,7-dihydroxy-6-(2-methoxypropan-2-yl)-8-(2-methylpropanoyl)-4-pentyl-2H-chromen-2-one

Chemical Properties


CAS No. :

Formula: C₂₂H₃₀O₆

Molecular Weight: 390.48

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

5,7-dihydroxy-6-(2-methoxypropan-2-yl)-8-(2-methylpropanoyl)-4-pentyl-2H-chromen-2-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.561 mL	12.8048 mL	25.6095 mL
5 mM	0.5122 mL	2.561 mL	5.1219 mL
10 mM	0.2561 mL	1.2805 mL	2.561 mL
50 mM	0.0512 mL	0.2561 mL	0.5122 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481