

## Compound TN12505

## Chemical Properties


CAS No. :

Formula: C<sub>18</sub>H<sub>20</sub>O<sub>10</sub>

Molecular Weight: 396.35

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description

4-(acetyloxy)-3-(((2E)-3-(3,4-dihydroxyphenyl)prop-2-enoyl)oxy)-1,5-dihydroxycyclohexane-1-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.523 mL	12.6151 mL	25.2302 mL
5 mM	0.5046 mL	2.523 mL	5.046 mL
10 mM	0.2523 mL	1.2615 mL	2.523 mL
50 mM	0.0505 mL	0.2523 mL	0.5046 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481