

## Compound TN12508

## Chemical Properties


CAS No. :

Formula: C47H74O18

Molecular Weight: 927.09

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(2S,3R,4R,5R,6R)-6-(((3R,6aR,6bS,8aS,12aS,14bR)-4,4,6a,6b,11,11,14b-heptamethyl-8a-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)carbonyl)-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,14,14a,14b-icosahydricen-3-yl]oxy)-4,5-dihydroxy-3-(((2S,3R,4S,5S)-3,4,5-trihydroxyoxan-2-yl]oxy)oxane-2-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0786 mL	5.3932 mL	10.7864 mL
5 mM	0.2157 mL	1.0786 mL	2.1573 mL
10 mM	0.1079 mL	0.5393 mL	1.0786 mL
50 mM	0.0216 mL	0.1079 mL	0.2157 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481