

Compound TN12511

Chemical Properties


CAS No. :

Formula: C₃₁H₂₈O₁₅

Molecular Weight: 640.55

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	[(2R,3R,4S,5R,6S)-4-(((2E)-3-(3,4-dihydroxyphenyl)prop-2-enoyl[oxy]-3,5-dihydroxy-6-(((2E)-3-(4-hydroxyphenyl)prop-2-enoyl[oxy]oxan-2-yl)methyl 3,4,5-trihydroxybenzoate)] is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5612 mL	7.8058 mL	15.6116 mL
5 mM	0.3122 mL	1.5612 mL	3.1223 mL
10 mM	0.1561 mL	0.7806 mL	1.5612 mL
50 mM	0.0312 mL	0.1561 mL	0.3122 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481