

Compound TN12528

Chemical Properties


CAS No. :

Formula: C37H52O5

Molecular Weight: 576.82

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(1S,2R,4aS,6aS,6bR,10R,11R,12aR)-11-(benzoyloxy)-10-hydroxy-1,2,6a,6b,9,9,12a-heptamethyl-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,12b,13,14b-icosahydricene-4a-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7336 mL	8.6682 mL	17.3364 mL
5 mM	0.3467 mL	1.7336 mL	3.4673 mL
10 mM	0.1734 mL	0.8668 mL	1.7336 mL
50 mM	0.0347 mL	0.1734 mL	0.3467 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481