

8-(hydroxymethyl)-4b,8-dimethyl-2-(propan-2-yl)-4b,5,6,7,8,8a,9,10-octahydrophenanthren-3-ol

Chemical Properties

CAS No. :


Formula: C₂₀H₃₀O₂

Molecular Weight: 302.46

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

**Biological Description**

Description

(4bS,8R)-8-(hydroxymethyl)-4b,8-dimethyl-2-(propan-2-yl)-4b,5,6,7,8,8a,9,10-octahydrophenanthren-3-ol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3062 mL	16.5311 mL	33.0622 mL
5 mM	0.6612 mL	3.3062 mL	6.6124 mL
10 mM	0.3306 mL	1.6531 mL	3.3062 mL
50 mM	0.0661 mL	0.3306 mL	0.6612 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481