

6-hydroxy-7,8-dimethoxy-3-phenyl-1H-isochromen-1-one

Chemical Properties


CAS No. :

Formula: C17H14O5

Molecular Weight: 298.29

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

6-hydroxy-7,8-dimethoxy-3-phenyl-1H-isochromen-1-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 3.3524 mL | 16.7622 mL | 33.5244 mL |
| 5 mM | 0.6705 mL | 3.3524 mL | 6.7049 mL |
| 10 mM | 0.3352 mL | 1.6762 mL | 3.3524 mL |
| 50 mM | 0.067 mL | 0.3352 mL | 0.6705 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481