

Compound TN12542

Chemical Properties


CAS No. :

Formula: C50H78O20

Molecular Weight: 999.15

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2S,3S,4S,5R,6R)-6-(((3S,4R,6aR,6bS,8aR,9R,11R,12aS,14bR)-9-(acetyloxy)-11-formyl-4-(hydroxymethyl)-4,6a,6b,8a,11,14b-hexamethyl-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,14,14a,14b-icosahydropicen-3-yl]oxy)-5-(((2S,3R,4S,5R,6S)-4,5-dihydroxy-6-(hydroxymethyl)-3-(((2S,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy)oxan-2-yl]oxy)-3,4-dihydroxyoxane-2-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0009 mL	5.0043 mL	10.0085 mL
5 mM	0.2002 mL	1.0009 mL	2.0017 mL
10 mM	0.1001 mL	0.5004 mL	1.0009 mL
50 mM	0.020 mL	0.1001 mL	0.2002 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481