

## Compound TN12555

## Chemical Properties


CAS No. :

Formula: C<sub>25</sub>H<sub>28</sub>O<sub>12</sub>

Molecular Weight: 520.49

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description

[(2R,3R,4S,5R,6S)-4-(acetyloxy)-6-{3,5-dihydroxy-2-[3-(4-hydroxyphenyl)propanoyl]phenoxy}-3,5-dihydroxyoxan-2-yl]methyl acetate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9213 mL	9.6063 mL	19.2127 mL
5 mM	0.3843 mL	1.9213 mL	3.8425 mL
10 mM	0.1921 mL	0.9606 mL	1.9213 mL
50 mM	0.0384 mL	0.1921 mL	0.3843 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481