

Compound TN12570

Chemical Properties


CAS No. :

Formula: C₄₀H₆₂O₁₅

Molecular Weight: 782.92

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	6-[(1R)-1-[(1S,3bS,7R,9S,9aR,11aS)-9-hydroxy-9a,11a-dimethyl-7-[[[(2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-((2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)methyl]oxan-2-yl]oxy]-1H,2H,3H,3aH,3bH,4H,6H,7H,8H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-1-yl]-1-hydroxyethyl]-3,4-dimethyl-5,6-dihydro-2H-pyran-2-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2773 mL	6.3863 mL	12.7727 mL
5 mM	0.2555 mL	1.2773 mL	2.5545 mL
10 mM	0.1277 mL	0.6386 mL	1.2773 mL
50 mM	0.0255 mL	0.1277 mL	0.2555 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481