

2-[[1,3-dihydroxy-1-(4-hydroxyphenyl)propan-2-yl]oxy]-6-(hydroxymethyl)oxane-3,4,5-triol

Chemical Properties


CAS No. :

Formula: C₁₅H₂₂O₉

Molecular Weight: 346.33

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2R,3R,4S,5S,6R)-2-[[1,3-dihydroxy-1-(4-hydroxyphenyl)propan-2-yl]oxy]-6-(hydroxymethyl)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8874 mL	14.4371 mL	28.8742 mL
5 mM	0.5775 mL	2.8874 mL	5.7748 mL
10 mM	0.2887 mL	1.4437 mL	2.8874 mL
50 mM	0.0577 mL	0.2887 mL	0.5775 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481