

Compound TN12594

## Chemical Properties


CAS No. :

Formula: C53H86O23

Molecular Weight: 1091.25

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(2S,3S,4S,5R,6R)-6-[[[(3S,4S,6aR,6bR,8aR,9R,11S,14bR)-9-hydroxy-4,11-bis(hydroxymethyl)-4,6a,6b,8a,11,14b-hexamethyl-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,14,14a,14b-icosahydricen-3-yl]oxy]-3,4-dihydroxy-5-[[[(2S,3R,4S,5S,6R)-5-hydroxy-6-(hydroxymethyl)-3-[[[(2S,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy]-4-[[[(2S,3R,4S,5S)-3,4,5-trihydroxyoxan-2-yl]oxy]oxan-2-yl]oxy]oxane-2-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9164 mL	4.5819 mL	9.1638 mL
5 mM	0.1833 mL	0.9164 mL	1.8328 mL
10 mM	0.0916 mL	0.4582 mL	0.9164 mL
50 mM	0.0183 mL	0.0916 mL	0.1833 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481