

## Compound TN12604

## Chemical Properties


CAS No. :

Formula: C<sub>40</sub>H<sub>60</sub>O<sub>16</sub>

Molecular Weight: 796.9

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	6-[1-(3a-hydroxy-9a,11a-dimethyl-9-oxo-7-(((2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(((2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)methyl)oxan-2-yl]oxy)-1H,2H,3H,3aH,3bH,4H,6H,7H,8H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-1-yl)-1-hydroxyethyl]-3,4-dimethyl-5,6-dihydro-2H-pyran-2-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2549 mL	6.2743 mL	12.5486 mL
5 mM	0.251 mL	1.2549 mL	2.5097 mL
10 mM	0.1255 mL	0.6274 mL	1.2549 mL
50 mM	0.0251 mL	0.1255 mL	0.251 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481