

Compound TN12611

Chemical Properties


CAS No. :

Formula: C₂₉H₃₈O₁₀

Molecular Weight: 546.61

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-hydroxy-2,3-dimethyl-4-((3,4,5,8-tetramethoxy-9,10-dimethyl-15,17-dioxatetracyclo[10.7.0.0 ^{2,7} .0 ^{14,18}]nonadeca-1(12),2(7),3,5,13,18-hexaen-19-yl)oxy)butanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8295 mL	9.1473 mL	18.2946 mL
5 mM	0.3659 mL	1.8295 mL	3.6589 mL
10 mM	0.1829 mL	0.9147 mL	1.8295 mL
50 mM	0.0366 mL	0.1829 mL	0.3659 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481