

Compound TN12619

Chemical Properties


CAS No. :

Formula: C₃₉H₆₄O₁₅

Molecular Weight: 772.93

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2R,3R,4S,5S,6R)-2-[[[(2S,4S,6R,7S,9S,13R,14R,16S)-6,14-dihydroxy-7,9,13-trimethyl-6-[(3R)-3-methyl-4-[[[(2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy}butyl]-5-oxapentacyclo[10.8.0.0 ^{2,9} .0 ^{4,8} .0 ^{13,18}]]icos-18-en-16-yl]oxy]-6-(hydroxymethyl)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2938 mL	6.4689 mL	12.9378 mL
5 mM	0.2588 mL	1.2938 mL	2.5876 mL
10 mM	0.1294 mL	0.6469 mL	1.2938 mL
50 mM	0.0259 mL	0.1294 mL	0.2588 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481