

## Compound TN12638

### Chemical Properties


CAS No. :

Formula: C<sub>27</sub>H<sub>44</sub>O<sub>7</sub>

Molecular Weight: 480.64

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description	(3a <sub>S</sub> ,5a <sub>R</sub> ,7 <sub>R</sub> ,8 <sub>S</sub> ,9a <sub>R</sub> ,11a <sub>R</sub> )-3a,7,8-trihydroxy-9a,11a-dimethyl-1-[(2 <sub>R</sub> ,3 <sub>R</sub> )-2,3,7-trihydroxy-6-methylheptan-2-yl]-1H,2H,3H,3aH,5H,5aH,6H,7H,8H,9H,9aH,9bH,10H,11H,11aH-cyclopenta[a]phenanthren-5-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0806 mL	10.4028 mL	20.8056 mL
5 mM	0.4161 mL	2.0806 mL	4.1611 mL
10 mM	0.2081 mL	1.0403 mL	2.0806 mL
50 mM	0.0416 mL	0.2081 mL	0.4161 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481